



# Knee Brace

The Knapp Knee Brace is more than a hinged knee brace....

It is so versatile and flexible, we call it a system.

Universal Sizing – NoSwet® Breathable Fabric  
Custom Shaped to Fit the Patient.

Hinges are offset in the center used in combination with pads, they will not irritate the knee condyle.

## Physician Inspired Bracing



### Unique Features

Constructed of Kuhl™ perforated neoprene for cool, comfortable compression

One brace should cover MCL & LCL sprains, meniscus injuries-even mild to moderate ACL & PCL Injuries!

2 Lengths effectively handle multiple instability diagnoses

3 Hinge styles increase customization

Hinges are fully covered for sports activity, yet easily accessible

Elasticized band in thigh strap ensures comfort

Ideal for pre-post-operative care of MCL/LCL and meniscal repairs

Longer length can be used for mild-moderate ACL/PCL sprains & strains and as an adjunct to patella fracture

Open popliteal eliminates bunching, improves range of motion and cools the back of the knee

Hinges are offset in the center to reduce condyle irritation.



### Natural Motion Restored

*“The primary function of a hinged knee brace is to increase the stability of the knee joint, yet allow flexion and extension.*

Dr. Thomas Knapp

**ZYNEX**  
MEDICAL™  
9990 Park Meadows Dr.  
Lone Tree, CO 80124  
Phone (866) 503-8921  
Fax (800) 495-6695  
www.zynexmed.com